


Product Info

| | | |
|-------------------------|--------------------------|---------------------------|
| Date: 17/3 - 2006 | Product info.: 4-2006 | Written by: Fredrik Delås |
| Item: New calf supports | Delivery: From this date | |

From this date we will start to deliver new calf supports.



The new calf support has gone through a "face lift" and is now made of moulded polyurethane. It is easy to clean, and has no plastic edges or cracks where dirt has a tendency to find its way. The first calf supports in sizes 43, 45 & 50 cm are going into production today. The new calf support will be implemented on all by the end of April.

| Item number | Item |
|-------------|---|
| 27964 | Calf support 35cm complete w/bracket |
| 27965 | Calf support 38-40cm complete w/bracket |
| 27966 | Calf support 43-45cm complete w/bracket |
| 27967 | Calf support 50cm complete w/bracket |
| 27968 | Calf support 55cm complete w/bracket |
| 27969 | Calf support 60cm complete w/bracket |