


Product Info

Date: 29/8 - 2006	Product info.: 8-2006	Written by: Fredrik Delås
Item: Foot plates	Delivery: From this date	

From this date we are supplying new foot plates for the Netti wheelchairs.



The new foot plate has been made more square to provide a larger contact area for the foot. The hinge has been changed to make the foot plate more stable. There is also added a spring under the foot plate to allow easier handling of the foot plate by the user. These foot plates will be mounted on all Netti chairs except for the Netti Mini models.

Item: Foot plates for Foot Support Universal with calf support	Item number
Foot plate 35cm with adjustment piece 245mm. With lock . Right	27850
Foot plate 38cm with adjustment piece 245mm. With lock. Right	27851
Foot plate 40cm with adjustment piece 245mm. With lock. Right	27852
Foot plate 43cm with adjustment piece 245mm. With lock. Right	27853
Foot plate 45cm with adjustment piece 245mm. With lock. Right	27854
Foot plate 50cm with adjustment piece 245mm. With lock. Right	27855
Foot plate 55cm with adjustment piece 245mm. With lock. Right	27856
Foot plate 60cm with adjustment piece 245mm. With lock. Right	27857
Foot plate 35cm with adjustment piece 245mm. With lock. Left.	27859
Foot plate 38cm with adjustment piece 245mm. With lock. Left.	27860
Foot plate 40cm with adjustment piece 245mm. With lock. Left.	27861
Foot plate 43cm with adjustment piece 245mm. With lock. Left.	27862
Foot plate 45cm with adjustment piece 245mm. With lock. Left.	27863
Foot plate 50cm with adjustment piece 245mm. With lock. Left.	27864
Foot plate 55cm with adjustment piece 245mm. With lock. Left.	27865
Foot plate 60cm with adjustment piece 245mm. With lock. Left.	27866

Item: Foot plates for standard Angle Adjustable Foot Fupports.	Item number*
Foot plate 35cm with adjustment piece 350mm. With lock. Right.	28621
Foot plate 38cm with adjustment piece 350mm. With lock. Right.	28622
Foot plate 40cm with adjustment piece 350mm. With lock. Right.	28623
Foot plate 43cm with adjustment piece 350mm. With lock. Right.	28624
Foot plate 45cm with adjustment piece 350mm. With lock. Right.	28625
Foot plate 50cm with adjustment piece 350mm. With lock. Right.	28626
Foot plate 55cm with adjustment piece 350mm. With lock. Right.	28627
Foot plate 60cm with adjustment piece 350mm. With lock. Right.	28628
Foot plate 35cm with adjustment piece 350mm. With lock. Left.	28649
Foot plate 38cm with adjustment piece 350mm. With lock. Right.	28650
Foot plate 40cm with adjustment piece 350mm. With lock. Right.	28651
Foot plate 43cm with adjustment piece 350mm. With lock. Right.	28652
Foot plate 45cm with adjustment piece 350mm. With lock. Right.	28653
Foot plate 50cm with adjustment piece 350mm. With lock. Right.	28654
Foot plate 55cm with adjustment piece 350mm. With lock. Right.	28655
Foot plate 60cm with adjustment piece 350mm. With lock. Right.	28656

Item: Foot plates for old Fixed Foot Supports & Foot Support Universal Standard	Item number**
Foot plate 35cm with adjustment piece 245mm. With heel strap. Right	27875
Foot plate 38cm with adjustment piece 245mm. With heel strap. Right	27876
Foot plate 40cm with adjustment piece 245mm. With heel strap. Right	27877
Foot plate 43cm with adjustment piece 245mm. With heel strap. Right	27878
Foot plate 45cm with adjustment piece 245mm. With heel strap. Right	27879
Foot plate 50cm with adjustment piece 245mm. With heel strap. Right	27880
Foot plate 55cm with adjustment piece 245mm. With heel strap. Right	27881
Foot plate 60cm with adjustment piece 245mm. With heel strap. Right	27882
Foot plate 35cm with adjustment piece 245mm. With heel strap. Left.	27899
Foot plate 38cm with adjustment piece 245mm. With heel strap. Left.	27900
Foot plate 40cm with adjustment piece 245mm. With heel strap. Left.	27901
Foot plate 43cm with adjustment piece 245mm. With heel strap. Left.	27902
Foot plate 45cm with adjustment piece 245mm. With heel strap. Left.	27903
Foot plate 50cm with adjustment piece 245mm. With heel strap. Left.	27904
Foot plate 55cm with adjustment piece 245mm. With heel strap. Left.	27905
Foot plate 60cm with adjustment piece 245mm. With heel strap. Left.	27906

Item: Foot plates for old Fixed Foot Supports & Foot Support Universal Standard	Item number***
Foot plate 35cm with adjustment piece 160mm. With heel strap. Right	27867
Foot plate 38cm with adjustment piece 160mm. With heel strap. Right	27868
Foot plate 40cm with adjustment piece 160mm. With heel strap. Right	27869
Foot plate 43cm with adjustment piece 160mm. With heel strap. Right	27870
Foot plate 45cm with adjustment piece 160mm. With heel strap. Right	27871
Foot plate 50cm with adjustment piece 160mm. With heel strap. Right	27872
Foot plate 55cm with adjustment piece 160mm. With heel strap. Right	27873
Foot plate 60cm with adjustment piece 160mm. With heel strap. Right	27874
Foot plate 35cm with adjustment piece 160mm. With heel strap. Left.	27891
Foot plate 38cm with adjustment piece 160mm. With heel strap. Left.	27892
Foot plate 40cm with adjustment piece 160mm. With heel strap. Left.	27893
Foot plate 43cm with adjustment piece 160mm. With heel strap. Left.	27894
Foot plate 45cm with adjustment piece 160mm. With heel strap. Left.	27895
Foot plate 50cm with adjustment piece 160mm. With heel strap. Left.	27896
Foot plate 55cm with adjustment piece 160mm. With heel strap. Left.	27897
Foot plate 60cm with adjustment piece 160mm. With heel strap. Left.	27898

Item: Foot plates for old Fixed Foot Supports & Foot Support Universal Standard	Item number****
Foot plate 35cm with adjustment piece 350mm. With heel strap. Right	27883
Foot plate 38cm with adjustment piece 350mm. With heel strap. Right	27884
Foot plate 40cm with adjustment piece 350mm. With heel strap. Right	27885
Foot plate 43cm with adjustment piece 350mm. With heel strap. Right	27886
Foot plate 45cm with adjustment piece 350mm. With heel strap. Right	27887
Foot plate 50cm with adjustment piece 350mm. With heel strap. Right	27888
Foot plate 55cm with adjustment piece 350mm. With heel strap. Right	27889
Foot plate 60cm with adjustment piece 350mm. With heel strap. Right	27890
Foot plate 35cm with adjustment piece 350mm. With heel strap. Left.	27907
Foot plate 38cm with adjustment piece 350mm. With heel strap. Left.	27908
Foot plate 40cm with adjustment piece 350mm. With heel strap. Left.	27909
Foot plate 43cm with adjustment piece 350mm. With heel strap. Left.	27910
Foot plate 45cm with adjustment piece 350mm. With heel strap. Left.	27911
Foot plate 50cm with adjustment piece 350mm. With heel strap. Left.	27912
Foot plate 55cm with adjustment piece 350mm. With heel strap. Left.	27913
Foot plate 60cm with adjustment piece 350mm. With heel strap. Left.	27914

- * Item numbers 28621-28628 & 28649-28656 replaces item numbers 22222-22227 & 22234-22259.
- ** Item numbers 27875-27882 & 27899-27906 replaces item numbers 22260-22273.
- *** Item numbers 27867-27874 & 27891-27898 replaces item numbers 22240-22255.
- **** Item numbers 27883-27890 & 27907-27914 replaces item numbers 22280-22295.