

## Product Info

Date: 26/9 - 2006	Product info.: 9-2006	Written by: Fredrik Delås
Item: New adjustment wheel	Delivery: From this date	

From this date we are starting to mount new adjustment wheels on foot support Universal on the following models:  
**Netti 4U Comfort CE & Netti 4U Comfort CED.**



Today we are mounting a bigger and more ergonomically shaped adjustment wheel on the Universal foot support. The colour is red to signal that it is to be adjusted.

This will in the first round be mounted on the foot support Universal used on the Netti 4U Comfort CE & Netti 4U Comfort CED models. In the near future these adjustment wheels will be used on all Universal foot supports in the Netti assortment. This will be implemented gradually during the next few months.

In centre of the adjustment wheel there is a black screw. This screw will on the first shipments be silver grey, and then we will go over to using a black screw.

**NB!** Remember to always turn the adjustment wheel all the way out before adjusting the angle.

Item number	Item	Model
29076	M8 Star wheel Ø60x30	NEW
23581	M8x45 Hex socket button head	NEW
26684	Star wheel M8 for foot & head support	OLD