

# PRODUCT INFORMATION

Date: 28/1 - 2008	Product info.: 1-2008	Written by: Fredrik Delås
Item: Oxy3 outer cover	Delivery: From this date	

## Oxy3 outer cover replaces Mars outer cover on seat cushion Sit.

To offer our customers an even more enhanced sitting comfort, we will now supply seat cushion Sit with outer cover Oxy3. The change starts from this date, but there will in a period of time be delivered cushions with Mars cover. This depends on numbers and sizes of cushions with Mars in stock.

Oxy3 is a textile made from 100% Polyester – with extraordinary good ventilating properties. Oxy3 has an open mesh structure, offering a high degree of ventilation, prevents heat build up and contributes to improvement of the micro climate in the cushion.

Oxy3 can be machine washed at 60°C and dried in a tumble drier at maximum 85°C.

The textile is flame retardant, durable and very elastic.

This textile is being used in aircrafts, within health care and in the fire sector and is therefore subdued to very strict usage- and documentation demands. The process of production is environmentally controlled and the textile can be recycled.

Vital Base and Alu Rehab is now working together towards that the Oxy3 should replace the Mars on all products in the Vital Base range of seating solutions.

In addition to the change of textile the Sit cushion has gone through a design improvement. The logo has been changed, and the core of the cushion is now rounded at the rear. This makes it easier to understand correct placement of the cushion in the wheelchair, and the rounded rear-end makes the cushion more versatile for usage in other wheelchairs. The core of the cushion is the same as before.

This change does NOT imply change of item numbers.



Sit Mars

Sit Oxy3