

Seat cushions from Alu Rehab

ALU REHAB™
Thinking ahead - moving forward

– Velcro on underside of seat cushions reintroduced



Following feedback and reports from the market, we have decided to reintroduce Velcro on the undersides of all seat cushions. We shall continue using friction cloth (info on friction cloth in Product Info no. 3/2008) as the base textile, but now in combination with Velcro. This combination will give improved stability and prevent the cushion sliding on the seat.

The change will have no effect on article number or price.

PRODUCT INFORMATION 2008/11

NO: 2008/11 | DELIVERY: FROM CURRENT DATE | www.alurehab.no

Date: 27.08.2008