

NEW PRODUCT:

Octopus drum brake wheels



NEW article numbers:

Drum brake wheel Octopus 22"

- art. no: 42078

Drum brake wheel Octopus 24"

- art. no: 42079

- Elegant new design
- Lighter weight
- Only supplier using such wheels

These wheels are 100% compatible with our existing model and there is no need to change any parts.

When compared to earlier models, these new drum brake wheels will make our chairs 275 grams lighter. Octopus has been tested and is just as strong as our current wheels.

Octopus drum brake wheels will replace our current model with effect from today's date, and they will be sent out as we gradually run out of stocks of the old model. They have NEW article numbers, but the price will remain the same.

PRODUCT INFORMATION 2009/1