

Accessories...

with inter-changeability!

A unique feature of the Netti wheelchairs is the vast amount of accessories available. Most of the accessories fit all the wheelchair models. We can also provide customized solutions according to individual needs!

HEAD SUPPORTS



Head support with sidesupport.



Head support small.



Head support large.

BACKREST CUSHIONS



Vital Base Smart.



Vital Base Stabil.




The backrest cushions are available in height 48 cm and 60 cm.

SEAT CUSHIONS



Vital Base Sit.



Vital Base Venti.



Vital Base Royal.

TRAYS



Universal trays.



Universal tray with padding.



Universal tray with lock and swingable.

WHEELS



Easy wheel for easier driving.



16" wheels with companion brake.



Different sizes of wheels with companion brake.

Why choose a Netti wheelchair?

Netti comfort wheelchairs are designed to give you stability and comfort. There is a wide choice of models within the Netti range and the design, combined with the inter-changeability of many parts, provides the possibility of adapting the wheelchair for your particular needs.

The manufacturer Alu Rehab has a policy of producing products which are functional, long lasting and cost effective for the users over the product life cycle.

We need to take care of our environment. Netti III comfort was the first wheelchair which has gone through EPD (Environment Product Declaration), which describes the impact a wheelchair has on the environment, from production to delivery to the customer.

All our wheelchairs are crash tested and approved according to ISO 7176/19.

For the future you can expect good and environmental friendly products from Alu Rehab.

GOOD COMFORT IS IN THE DETAILS!



Comfort wheelchairs

– designed to give you stability and comfort



PRODUCER:

ALU REHAB™

Thinking ahead - moving forward

ALU REHAB™

Thinking ahead - moving forward

Comfort wheelchairs

Designed to give you stability and comfort

Netti  **III comfort**
– the complete solution for complex needs!



Technical specifications:

Item nr.	Seat width	Seat depth (min.–max.)	Seat height (min.–max.)
20685	35 cm	40–50 cm	40–48 cm
20690	38 cm	40–50 cm	40–48 cm
20686	40 cm	40–50 cm	40–48 cm
20689	43 cm	40–50 cm	40–48 cm
20687	45 cm	40–50 cm	40–48 cm
20688	50 cm	40–50 cm	40–48 cm

Max user weight: 145 kg

Netti  **III comfort HD**
– the complete solution – in a larger size!



Technical specifications:

Item nr.	Seat width	Seat depth (min.–max.)	Seat height (min.–max.)
25369	50 cm	40–50 cm	40–48 cm
20691	55 cm	40–50 cm	40–48 cm
20694	60 cm	40–50 cm	40–48 cm

Max user weight: 160 kg

Netti  **comfort | CED**
– makes everyday life easier – with extra features!



Technical specifications:

Item nr.	Seat width	Seat depth (min.–max.)	Seat height (min.–max.)
28423	35 cm	42–50 cm	46–50 cm
28424	40 cm	42–50 cm	46–50 cm
40608	43 cm	42–50 cm	46–50 cm
28370	45 cm	42–50 cm	46–50 cm
28425	50 cm	42–50 cm	46–50 cm
40616	55 cm	42–50 cm	46–50 cm
40617	60 cm	42–50 cm	46–50 cm

Max user weight: 140 kg

Netti  **comfort | CEDS**
– the little chair that makes everyday life easier!



Technical specifications:

Item nr.	Seat width	Seat depth (min.–max.)	Seat height (min.–max.)
40875	35 cm	37–45 cm	44–48 cm
40876	40 cm	37–45 cm	44–48 cm
40877	43 cm	37–45 cm	44–48 cm
40878	45 cm	37–45 cm	44–48 cm
40879	50 cm	37–45 cm	44–48 cm

Max user weight: 130 kg

Netti  **comfort | FA**
– being totally adjustable – means less is more!



Technical specifications:

Item nr.	Seat width	Seat depth (min.–max.)	Seat height (min.–max.)
42076	35–44 cm	39–51 cm	40–49 cm
42077	44–53 cm	39–51 cm	40–49 cm
41300	35–44 cm	39–51 cm	40–49 cm
with drumbrake			
41301	44–53 cm	39–51 cm	40–49 cm
with drumbrake			

Max user weight: 140 kg

Netti  **comfort | CE**
– the simple solution!



Technical specifications:

Seat width	Seat depth (min.–max.)	Seat height (min.–max.)
35 cm	42–50 cm	46–50 cm
40 cm	42–50 cm	46–50 cm
45 cm	42–50 cm	46–50 cm
50 cm	42–50 cm	46–50 cm

Max user weight: 125 kg. Contact us for item nr.

Netti  **comfort | CES**
– the simple solution in a smaller size!



Technical specifications:

Seat width	Seat depth (min.–max.)	Seat height (min.–max.)
35 cm	37–45 cm	44–48 cm
40 cm	37–45 cm	44–48 cm
45 cm	37–45 cm	44–48 cm
50 cm	37–45 cm	44–48 cm

Max user weight: 125 kg. Contact us for item nr.

Netti  **I comfort**
– for self-propelling users!



Technical specifications:

Item nr.	Seat width	Seat depth (min.–max.)	Seat height (min.–max.)
24040	35 cm	40–50 cm	32–48 cm
24041	38 cm	40–50 cm	32–48 cm
24042	40 cm	40–50 cm	32–48 cm
24043	43 cm	40–50 cm	32–48 cm
24044	45 cm	40–50 cm	32–48 cm
24045	50 cm	40–50 cm	32–48 cm
24046	55 cm	40–50 cm	32–48 cm
24047	60 cm	40–50 cm	32–48 cm

Max user weight: 140 kg

Netti  **svipp comfort**
– for people on the move!



Technical specifications:

Item nr.	Seat width	Seat depth (min.–max.)	Seat height (min.–max.)
24160RC	35 cm	40–50 cm	40–48 cm
24161RC	38 cm	40–50 cm	40–48 cm
24162RC	40 cm	40–50 cm	40–48 cm
24163RC	43 cm	40–50 cm	40–48 cm
24164RC	45 cm	40–50 cm	40–48 cm
24165RC	50 cm	40–50 cm	40–48 cm
24166RC	55 cm	40–50 cm	40–48 cm

Max user weight: 125 kg

Netti  **mini**
– grows with the child!



Technical specifications:

Item nr.	Seat width	Seat depth (min.–max.)	Seat height (min.–max.)
23790	25 cm	32–38 cm	42–45 cm
23791	30 cm	32–38 cm	42–45 cm
23792	35 cm	32–38 cm	42–45 cm

Max user weight: 90 kg